

### Guinea pig polyclonal antibody packages

Customer will get:	Peptide (synthesis)	Antiserum	Preimmune serum	ELISA testing and results	IgG-specific purification	Epitope-specific purification
Package A	5-10 mg of crude up to 18 aa long peptide, MS and HPLC data	~7 ml per g. pig (from 1 of 2 g. pigs)	~ 0.1 ml per g. pig	+	-	From 5 ml of antiserum (standard)
Package B	5-10 mg of crude up to 18 aa long peptide, MS and HPLC data	~7 ml per g. pig (from 1 of 2 g. pigs)	~ 0.1 ml per g. pig	+	5 mg of pure antibodies	-
Package C	5-10 mg of crude up to 18 aa long peptide, MS and HPLC data	~7 ml per g. pig	~ 0.1 ml per g. pig	+	-	-
Package D	-	~7 ml per g. pig (from 1 of 2 g. pigs)	~ 0.1 ml per g. pig	+	-	From 5 ml of antiserum (standard)
Package E	-	~7 ml per g. pig (from 1 of 2 g. pigs)	~ 0.1 ml per g. pig	+	5 mg of pure antibodies	-
Package F	-	~7 ml per g. pig	~ 0.1 ml per g. pig	+	-	-